



EasyVac Vacuum Pumps



The Ai EasyVac sliding vane rotary vacuum pumps are our most economical option. These pumps utilize high grade lubricating oil and come with an exhaust filter, making it quieter to run. Available sizes: 1.8 CFM single stage, 7 CFM single stage, and 9 CFM dual stage.

Model	EasyVac 1.8	EasyVac 7	EasyVac 9
Electrical Requirements	110V 60Hz or 220V 50/60Hz single phase ¼ HP, 190W	110V 60Hz or 220V 50/60Hz single phase ½ HP, 370W	110V 60Hz or 220V 50/60Hz single phase 750W (1H)
Pumping Rate	2 cfm (110V 60Hz), 1.8 cfm (220V 50/60Hz)	7 cfm (110V 60Hz), 6 cfm (220V 50/60Hz)	9 cfm (110V 60Hz), 8 cfm (220V 50/60Hz)
Vacuum Level	150 micron/milliorr	150 micron/milliorr	15 micron/milliorr
Exhaust Connection	KF25/NW25 flange		
Vacuum Connection	3/8" hose barb		KF25/NW25 flange
Vacuum Hose Type	3/8" white premium food grade silicone vacuum hose		KF25 stainless steel vacuum bellow
Rotary Vanes	single stage		Dual stage
Oil Volume	250 ml / 0.26 quart	415 ml / 0.44 quart	560 ml / 0.6 quart
Motor Rate	1720 rpm		
Working Temperature	-5 to 60° C		